

1 Full-Length Paper (Revised)

2

3 Title: Evaluation of a Loop-Mediated Isothermal Amplification (LAMP) Method as a

4 Diagnostic Tool of Zoonotic Simian malaria parasite *Plasmodium knowlesi* Infection

5

6 Running title: A LAMP METHOD FOR PLASMODIUM KNOWLESI INFECTION

7 Authors: Hiroshi Iseki,<sup>1</sup> Satoru Kawai,<sup>2\*</sup> Nobuyuki Takahashi<sup>1</sup>,<sup>1</sup> Makoto Hirai,<sup>3</sup>

8 Kazuyuki Tanabe,<sup>4</sup> Naoaki Yokoyama,<sup>1</sup> Ikuo Igarashi<sup>1</sup>

9

10 *National Research Center for Protozoan Diseases, Obihiro University of*  
11 *Agriculture and Veterinary Medicine, Obihiro, Hokkaido, Japan<sup>1</sup>; Center for*  
12 *Tropical Medicine and Parasitology, Dokkyo Medical University, Mibu,*  
13 *Tochigi, Japan<sup>2</sup>; Department of Parasitology, Gunma University Graduate*  
14 *School of Medicine, Maebashi, Gunma, Japan<sup>3</sup>; Laboratory of Malariology,*  
15 *Research Institute for Microbial Diseases, Osaka University, Suita, Osaka,*  
16 *Japan<sup>4</sup>*

17

18 \*Corresponding author. Mailing address: Satoru Kawai

19 Center for Tropical Medicine and Parasitology,  
20 Dokkyo Medical University,  
21 Mibu, Tochigi, 321-0293 Japan.  
22 Phone: 81 282 87 2134. Fax: 81 282 86 6431.  
23 Email:skawai@dokkyomed.ac.jp

## ABSTRACT

Loop-mediated isothermal amplification (LAMP) is a novel method that rapidly amplifies target DNA with high specificity under isothermal conditions. It has been applied as a diagnostic tool for several infectious diseases including viral, bacterial, and parasitic diseases. In the present study, we developed a LAMP method for the molecular diagnosis of *Plasmodium knowlesi* infection (*PkLAMP*) and evaluated its sensitivity, specificity, and clinical applicability. We designed three sets of *PkLAMP* primers for the species-specific  $\beta$ -tubulin gene. The primer sets for *PkLAMP* specifically amplified the autologous DNA extracts of *P. knowlesi*, and the sensitivity of the test was 100-fold that of single-PCR assay. These results indicate that our *PkLAMP* method can be used to efficiently distinguish between *P. knowlesi* and other malaria parasites. To evaluate the feasibility of using *in vivo* materials, comparisons of *PkLAMP* and the conventional nested PCR (nPCR) method and microscopic examination were made with blood samples from two experimentally infected monkeys. These studies showed that *PkLAMP* can be identified in the infectious course of *P. knowlesi* much earlier than with nPCR and microscopy. Moreover, the detection performance of *PkLAMP* using whole blood as the template was identical to that of *PkLAMP* when genomic DNA extracts were used. These results suggest that the *PkLAMP* method is a promising tool of molecular diagnosis of *P. knowlesi* infection in endemic areas.

## 1    **Introduction**

2            Naturally acquired human infections with a macaque malaria parasite, *Plasmodium*  
3    *knowlesi*, have now been referred to as the fifth human malaria (4, 17). In fact, recent studies  
4    have shown that naturally occurring *P. knowlesi* malaria cases are not rare and are widely  
5    distributed in Southeast Asia, particularly in forested areas inhabited by the natural macaque  
6    host and vectors such as the *Anopheles leucophyrus* group (4, 5, 16).

7            Until recently, numerous cases of *P. knowlesi* infections in humans may have been  
8    misdiagnosed as ordinary *P. malariae* malaria (4, 5, 16), since the morphological  
9    characteristics of the blood stages of *P. knowlesi* parasites are similar to those of *P. malariae*,  
10   and it can be easily misidentified as *P. malariae* on microscopic examination (16). Moreover,  
11   our recent study showed that some commercial rapid malaria diagnostic tests based on the  
12   detection of parasite lactate dehydrogenase enzyme (pLDH) are unable to distinguish between  
13   human malaria parasites and *P. knowlesi* since certain antibodies to pLDH that were thought to  
14   be specific for *P. falciparum* and *P. vivax* also bind to *P. knowlesi* (9). Although the  
15   development of a PCR diagnostic method has been essential to solving these problems of  
16   misdiagnosis, PCR assays are not a simple method of detection and are not a viable option for  
17   routine diagnosis.

18            Loop-mediated isothermal amplification (LAMP) has been developed as a novel method

1 to amplify DNA with high specificity and simplicity (13). It consists simply of incubating a  
2 mixture of the target gene, four or six different primers, *Bst* DNA polymerase, and substrates.  
3 The significant advantages of the LAMP method are (i) high amplification efficiency under  
4 isothermal conditions (63 to 65°C) and (ii) visual judgment based on the turbidity or  
5 fluorescence of the reaction mixture, which is kept in the reaction tube (10, 12). LAMP has  
6 thus emerged as a powerful tool to facilitate genetic testing for the rapid diagnosis of several  
7 infectious diseases including viral, bacterial, and parasitic diseases (8, 11). Although the  
8 detection performances of LAMP for four human malaria parasites have been assessed in  
9 clinical and epidemiological settings, the LAMP method has not yet been evaluated for the  
10 diagnosis of *P. knowlesi* infection (3, 7, 14). In the present study, we developed a LAMP  
11 method for diagnosis of *P. knowlesi* infection (*PkLAMP*) and evaluated its sensitivity,  
12 specificity, and clinical applicability using blood samples obtained from experimentally *P.*  
13 *knowlesi*-infected monkeys.

14

## 15 **Materials and Methods**

### 16 **Specific primers of *PkLAMP***

17 The LAMP method requires a set of four specific primers: a forward inner primer (FIP),  
18 a backward inner primer (BIP), and two outer primers (F3 and B3), which recognize a total of

1 six distinct nucleotide sequences (B1, B2, B3, F1, F2, and F3) on the target gene (10, 12, 13).  
2 Since it has been demonstrated that additional loop primers increase the amplification  
3 efficiency, loop primers for each target gene were also synthesized. The specific primers for *P.*  
4 *knowlesi* were designed against species-specific  $\beta$ -tubulin gene sequences (GenBank  
5 accession number: AY639984) (Fig. 1A). For easy confirmation of the amplified sequences,  
6 we modified FIP and BIP by inserting a restriction enzyme (*Eco* RI) cleavage site between the  
7 F1 complementary and F2 and between the B1 complementary and B2, respectively, as shown  
8 in Figure 1B.

#### 10 ***Pk*LAMP procedures**

11 The *Pk*LAMP reaction was performed as described previously (10, 12, 13). Briefly, the  
12 reaction was performed in 25  $\mu$ l of a mixture containing 1  $\mu$ l of the extracted DNA template,  
13 40 pmol each of the FIP and BIP primers, 5 pmol each of the F3 and B3 primers, 20 pmol each  
14 of the forward loop primer FLP and backward loop primer BLP, and 1  $\mu$ l of Fluorescent  
15 Detection Reagent (Eiken Chemical Co., Ltd. Tokyo, Japan) with LoopAmp DNA  
16 Amplification kit (Eiken Chemical Co., Ltd. Tokyo, Japan). The *Pk*LAMP reaction was  
17 performed as described above with each of the specific primers. In a conventional heat block,  
18 the mixture was incubated at 66 °C (47-72 °C was also tested) for 60 min and the reaction was

1 then terminated by heating the mixture at 80 °C for 5 min for termination. For the initial  
2 validation study, *PkLAMP* was confirmed with real-time monitoring of increase of turbidity  
3 using a Loopamp real-time turbidimeter (LA-200; Teramecs, Kyoto, Japan). To confirm the  
4 amplified DNA products of each parasite, 1 µg/µl of the product was digested with the *Eco* RI  
5 at 37 °C for 1 h. The non-treated and *Eco* RI digested LAMP product were subjected to  
6 electrophoresis on a 2% agarose gel and then visualized under ultraviolet (UV) light after  
7 staining with ethidium bromide (SIGMA). Digested LAMP DNA products were purified after  
8 2% agarose gel electrophoresis and then cloned into a pCRII cloning vector using a TA  
9 Cloning Kit (Invitrogen, Carlsbad, CA, USA). The nucleotide sequences of inserts were  
10 determined using a Big Dye Terminator Kit (Applied Biosystems Japan, Ltd.) with an  
11 automated DNA sequencer (ABI PRISM 3100 genetic analyzer, Applied Biosystems Japan,  
12 Ltd.). The Genetyx 7 package (Software Development Co., Ltd., Tokyo, Japan) was used to  
13 align the determined sequences. For the challenge infections, the amplified products in the  
14 reaction tube were directly detected with the naked eye using Loopamp fluorescent detection  
15 reagent (Eiken Chemical Co., Ltd.) according to the manufacturer's instructions.

16

#### 17 **Specificity of *PkLAMP* primers**

18 Specificity of the *PkLAMP* primers was tested using genomic DNA (gDNA) of various

1 *Plasmodium* species in a gel electrophoresis and fluorescent analysis. The gDNAs of *P.*  
2 *falciparum*, *P. vivax*, *P. malariae*, and *P. ovale* were kindly provided by Dr. Takefumi Tsuboi  
3 of Ehime University of Japan. The blood samples infected with *P. inui*, *P. simiovale*, *P. fieldi*,  
4 *P. fragile*, *P. hylobati*, and *P. gonderi* were obtained from American Type Culture Collection  
5 (ATCC) and gDNAs of these parasites were extracted from frozen infected blood by a  
6 QIAamp DNA Blood Mini Kit (QIAGEN, Tokyo, Japan), according to the manufacturer's  
7 instruction. The *P. coatneyi*- and *P. cynomolgi*-infected blood samples were obtained from  
8 experimentally infected monkeys, and were subjected to DNA extraction with the QIAamp  
9 DNA Blood Mini Kit. These purified DNA samples were used as templates for the subsequent  
10 *PkLAMP* and single-PCR assays. As a negative control, DNA extracted from normal monkey  
11 blood was prepared as described above.

12

### 13 **Sensitivity tests for *PkLAMP* and single-PCR**

14 For sensitivity tests, the *PkLAMP* reaction was tested using 10-fold serial dilutions of  
15 plasmid DNA containing the target sequence by cloning from *P. knowlesi* H strain genomic  
16 DNA and compared against results of the single-PCR assay using F3 and B3 primers. PCR  
17 amplification was performed in 25  $\mu$ l of a mixture containing 1  $\mu$ l of the extracted DNA  
18 template, 50 pmol of each primer, 200  $\mu$ M of each dNTP, and 1.25 U of Taq Gold DNA

1 polymerase (Applied Biosystems, Foster City, CA, USA) in a PCR buffer (Applied  
2 Biosystems). The reaction was performed at 35 cycles under the following conditions: 10 min  
3 at 95 °C to activate the Taq Gold DNA polymerase, 1 min of denaturation at 94 °C, 1 min of  
4 annealing at 60 °C, 1 min of extension at 72 °C, and 10 min of final extension at 72 °C in a  
5 Gene Amp PCR system 9700 (Applied Biosystems). The PCR products were subjected to  
6 agarose gel electrophoresis and then visualized as described above.

7

#### 8 **Evaluation of *Pk*LAMP using blood samples from infected monkeys**

9 *Pk*LAMP was practically evaluated for fluorescence detection of *P. knowlesi* target DNA  
10 using blood samples obtained from experimental *P. knowlesi*-infected monkeys. Two monkeys,  
11 J58 (male) and J64 (male), which were three-year-old Japanese macaques (*Macaca fuscata*)  
12 weighing 4.2 kg and 4.7 kg, respectively, were used in this experiment. Both monkeys were  
13 second-generation offspring bred in captivity. The investigators adhered to the Guidelines for  
14 the Use of Experimental Animals authorized by the Japanese Association for Laboratory  
15 Animal Science. Monkey J58 was inoculated intravenously with  $1 \times 10^8$  fresh *P. knowlesi* H  
16 strain (ATCC No. 30158) parasitized red blood cells (PRBCs) obtained from another infected  
17 Japanese macaque. Monkey J64 was inoculated intravenously with frozen *P. knowlesi* Hackeri  
18 strain (ATCC No. 30153) infected blood obtained from ATCC. After infection, Giemsa-stained



1 thin blood films were prepared daily from peripheral blood obtained by earprick, and  
2 parasitemia in the infected monkeys was monitored by microscopic examination. Heparinized  
3 blood samples for *Pk*LAMP assay were obtained daily from the infected monkeys during the  
4 course of infection. The infected blood samples were subjected to DNA extraction with a  
5 QIAamp DNA Blood Mini Kit (QIAGEN) as described above. The DNA extracts and whole  
6 blood samples were frozen at  $-80^{\circ}\text{C}$  until use.

#### 7 8 **Comparison of *Pk*LAMP and nested PCR using DNA extracts and whole blood as** 9 **template**

10 We compared the sensitivities of *Pk*LAMP and conventional nested PCR (nPCR) assays  
11 using DNA extract of *P. knowlesi* and whole blood obtained from two infected monkeys  
12 during the course of infection. The nPCR assay, based on the *Plasmodium* DNA sequence of  
13 the small subunit ribosomal RNA (SSUrRNA) genes, was performed according to a standard  
14 protocol as described previously (15). Nest 1 reaction was carried out in a 50 $\mu\text{l}$  reaction  
15 mixture containing 2x PCR master mix (Ampli Taq Gold PCR Master Mix, Applied  
16 Biosystems, New Jersey USA), 250nM of each primer (rPLU1 and rPLU5) (15) and 2  $\mu\text{l}$  of  
17 DNA template. The reaction mixture for Nest 1 PCR amplification was placed in a thermal  
18 cycler (TP600, Takara Bio inc., Shiga, Japan) at  $95^{\circ}\text{C}$  for 5 min for initial denaturation. This

1 was followed by 40 cycles of 94 °C for 30 s and 55 °C for 60 s and 72 °C for 120 s for  
2 amplification, and then 72 °C for 10 min for final extension. Nest 2 PCR amplification was  
3 performed in 20µl reaction mixture containing 2x PCR master mix (Applied Biosystems),  
4 250nM of each primer (Pmk8 and Pmkr9) (16) and 2 µl of the nest 1 PCR products used as  
5 DNA templates. A reaction mixture for Nest 2 PCR amplification was placed in a thermal  
6 cycler (TP600) at 95 °C for 5 min for initial denaturation. This was followed by 40 cycles of  
7 94 °C for 30 s and 60 °C for 60 s and 72 °C for 60 s for amplification, and then 72 °C for 10  
8 min for final extension. Nest 2 PCR products were electrophoresed separately on 2% agarose  
9 gel and illuminated with UV light.

10

## 11 **Results**

### 12 **Specificity of *PkLAMP* primers**

13 The specificity of *PkLAMP* primers was investigated by using various *Plasmodium*  
14 gDNAs as template for *PkLAMP*. As shown in Figure 2A, a typical ladder pattern was  
15 detected in *P. knowlesi* (lane 1) but not from the DNAs of other *Plasmodium* species (Fig. 2A).  
16 Moreover, fluorescent detection was also specifically obtained in the reaction tube including  
17 gDNA of *P. knowlesi* as shown in Figure 2B. The sizes of *PkLAMP* fragments digested by  
18 *Eco* RI were identical with the predicted sizes for the parasite (data not shown). To evaluate

1 the accuracy and robustness of the LAMP method, the *Pk*LAMP reaction was carried out in a  
2 water bath at 47-72 °C separately. Positive ladder patterns were observed at 48-71 °C, and  
3 strongly at 56-70 °C. These findings demonstrated that a set of species-specific primers was  
4 highly specific for the detection of the corresponding parasite in *Pk*LAMP. To confirm the  
5 nucleotide sequences of LAMP products, the amplified and digested DNA products were  
6 purified from the positive controls and cloned into a vector. The determined sequences of the  
7 DNA fragments were completely identical to the reported ones (data not shown) (*P. knowlesi*,  
8 AY639984).

#### 10 **Sensitivity of *Pk*LAMP reaction**

11 To examine the sensitivity of *Pk*LAMP, three detection methods of *Pk*LAMP were  
12 compared with conventional single-PCR using two outer primers; F3 and B3 for the detection  
13 of *P. knowlesi*  $\beta$ -tubulin gene. As shown in Figure 3A, amplification by real-time *Pk*LAMP  
14 was obtained in reaction tubes containing from  $10^8$  to  $10^2$  copies/ $\mu$ l of the DNA template for a  
15 60-min reaction with a turbidity assay. On gel electrophoresis analysis, the amplified products  
16 also showed ladder-like patterns from  $10^8$  to  $10^2$  copies/ $\mu$ l (Fig. 3B). The amplified products  
17 in these positive reaction tubes were also visually detectable using Loopamp fluorescent  
18 detection reagent, as shown in Figure 3C. In contrast, the limit of detection for PCR using F3

1 and B3 primers was  $10^8$  to  $10^4$  copies/ $\mu$ l (Fig. 3D). Therefore, it appeared that the sensitivity  
2 of the *Pk*LAMP, regardless of the detection methods, was 100-fold higher than that of  
3 single-PCR assay.

#### 4 5 **Evaluation of *Pk*LAMP and nPCR using DNA extracts and whole blood samples as a** 6 **template**

7 The course of *Macaca* monkeys experimental infected with *P. knowlesi* was monitored by  
8 *Pk*LAMP and nPCR reaction for the detecting parasite DNA (Table 1). Both monkeys infected  
9 with *P. knowlesi* developed a fulminating acute infection, and finally became lethargic and  
10 severely withdrawn just before autopsy. In monkey J58 inoculated with fresh PRBCs of *P.*  
11 *knowlesi* H strain, the parasites in the peripheral blood were first detected by microscopy on  
12 day 1; parasite densities then increased to around 10% within 3 days after infection. *P.*  
13 *knowlesi* DNA could be detected by *Pk*LAMP as well as nPCR assay on all days during the  
14 course of infection (Table 1). In monkey J64 inoculated with frozen PRBCs of *P. knowlesi*  
15 Hackeri strain, the parasites were first detected by microscopy on day 6; parasite densities then  
16 increased sharply to around 58% within 9 days after infection. *P. knowlesi* DNA could be  
17 detected by *Pk*LAMP throughout the course of infection, while the earliest detection of  
18 parasite DNA by the nPCR assay was the 3rd day after infection (Table 1).

1        We also compared the amplification efficiency of *PkLAMP* and nPCR using frozen whole  
2        blood as template. As shown in Table 1, *PkLAMP* could amplify the target from whole blood  
3        with similar efficiency to DNA extracts throughout the course of infection. These results  
4        clearly indicate that *PkLAMP* could detect even the target DNA from non-purified whole  
5        blood. In contrast, nPCR assay using whole blood from J58 and J64 could amplify parasite  
6        DNA only on day 3 and day 9, respectively, when parasite densities were markedly increased  
7        in the blood (Table 1).

8

## 9        **Discussion**

10        The diagnosis of malaria at regional clinics in endemic areas has mainly been performed  
11        by microscopic examination of blood smears because of its ease and rapid application.  
12        However, the morphology of the asexual stages of the zoonotic simian *Plasmodium* parasites  
13        substantially resembles to that of human parasites, particularly on thick blood-films, and  
14        laboratory technicians are trained to recognize only the four species of human parasites (16).  
15        In fact, numerous human cases of *P. knowlesi* infection have been misdiagnosed by  
16        microscopy as *P. malariae* due to their morphological similarities (4, 5, 16). The application of  
17        DNA amplification to the diagnosis of malaria can solve these problems. Amplification of  
18        parasite DNA using a specific PCR has been applied to various *Plasmodium* species including

1 four human malarial parasites and *P. knowlesi* (3, 7, 14). However, despite the excellent  
2 specificity and sensitivity of PCR and real-time PCR, these methods require complicated  
3 procedures and sophisticated instrumentation such as a thermal cycler, and they are often  
4 impracticable under conditions requiring field diagnosis. In this regard, the LAMP method has  
5 the advantages of simplicity, specificity, and sensitivity of reaction compared to other  
6 molecular diagnostic methods. It is thus the LAMP method is a promising candidate for wide  
7 use in regional clinics and under field conditions.

8 In the present study, we presented the successful development of a LAMP method for  
9 detecting *P. knowlesi* infection, using a primer set that targets the  $\beta$ -tubulin gene of parasites.  
10 The specificity of the primers was evaluated using nine species of simian malaria parasites and  
11 four species of human malaria parasites. The results showed that the primer set for *PkLAMP*  
12 amplified only the autologous DNA samples of *P. knowlesi* in typical ladder bands. In contrast,  
13 no ladder bands were obtained from any other control. These findings indicate that this primer  
14 set is specific for *P. knowlesi* and can be used to examine for knowlesi malaria as well as  
15 distinguish between it and other types of malaria. The sensitivity of the test was evaluated, and  
16 results showed that *PkLAMP* was 100-fold more sensitive than single-PCR assay using F3 and  
17 B3 primers. Moreover, the present study showed that an isothermal reaction time of 1hr was  
18 enough to amplify  $10^9$  of the target DNA in reaction tubes containing from  $10^8$  to  $10^2$  copies/ $\mu$ l

1 of the DNA template, and that results could be easily judged by visual inspection of the  
2 turbidity or fluorescence of the reaction mixture (10, 13). These results suggest that the  
3 *PkLAMP* assay is reliable and useful for the diagnosis of *knowlesi* malaria.

4 To evaluate the feasibility of using *in vivo* materials, comparisons of *PkLAMP* and the  
5 conventional nested PCR method and microscopic examination were made with blood samples  
6 from two infected monkeys. These studies validated *PkLAMP* as an alternative molecular  
7 diagnostic tool, which can be used in the diagnosis of early and advanced infections of *P.*  
8 *knowlesi*. Early species identification in the diagnosis for malaria is very important in  
9 preventing disease progression. In particular, early identification of *P. knowlesi* infection is  
10 essential, since the unique 24 hr asexual replication cycle among human and simian malaria  
11 parasites can rapidly result in high levels of parasitemia with a fatal outcome in humans (4, 5).  
12 Although nPCR and sequencing have been applied to species identification for malaria  
13 diagnosis, a more rapid diagnostic test such as *PkLAMP* would be a convenient and powerful  
14 tool for enabling the delivery of prompt and adequate medical treatment.

15 The present study also assessed the detection performance of *PkLAMP* with different  
16 DNA template preparations including frozen whole blood or genomic DNA extracts. The  
17 detection efficiency of *PkLAMP* using whole blood was identical to that of *PkLAMP* when  
18 gDNA extracts were used as the template. However, the detection performance of nPCR using

1 the whole blood templates was quite poor. It appears that this is due to blood components such  
2 as myoglobin, hem-blood protein complexes, and immunoglobulin G that inactivate *Taq* DNA  
3 polymerase used in standard PCR (2). In contrast, such inhibitors do not affect the *Bst*  
4 polymerase used in LAMP (6). According to previous reports, the specificity and sensitivity of  
5 detection appear to be unaffected by LAMP processing conditions or sample type, including  
6 whole blood, filter paper or card-processed blood, serum, sputum, and crudely processed tissue  
7 samples (8). Furthermore, Poon et al. have reported that *P. falciparum* DNA was detected by  
8 LAMP using a promising simple DNA template preparation method from heat-treated blood  
9 (14). Further improvement of template production methods for *Pk*LAMP will be required to  
10 optimize and simplify template preparation.

11 In conclusion, *Pk*LAMP can be considered as an efficient candidate for the molecular  
12 diagnosis of *P. knowlesi* infection in endemic areas. Thekisoe et al. reported that LAMP  
13 reagents are stable at ambient temperature for up to 2 weeks (17). In addition, a recent study of  
14 the LAMP method showed that it is able to detect both *Plasmodium* oocysts and sporozoites  
15 from an “all-in-one” template using whole mosquito bodies (1). These observations further  
16 emphasize the potential usefulness of the LAMP method as a diagnostic and new  
17 epidemiological surveillance tool for malaria. Our studies will also provide a powerful method  
18 for the diagnosis and monitoring of *P. knowlesi* infection in the field.



## 1    **Acknowledgments**

2        We are grateful to Mayumi Ohshita and Mayu Tanaka at Center for Tropical Medicine  
3    and Parasitology, Dokkyo Medical University for technical support. This study was supported  
4    in part by the following three grants; (i) A grant from the program for Promotion of  
5    Fundamental Studies in Health Sciences (No. 04-09) from the National Institute of Biomedical  
6    Innovation, (ii) A grant from the Promotion of Basic Research Activities for Innovative  
7    Biosciences program (PROBRAIN), and (iii) A grant from the 21<sup>st</sup> Century COE Program  
8    (A-1) of the Ministry of Education, Culture, Sports, Science, and Technology of Japan.

1   **References**

- 2   1.   **Aonuma H, Suzuki M, Iseki H, Perera N, Nelson B, Igarashi I, Yagi T, Kanuka H,**  
3       **Fukumoto S.** 2008. Rapid identification of *Plasmodium*-carrying mosquitoes using  
4       loop-mediated isothermal amplification. *Biochem Biophys Res Commun.* **376**: 671-675.
- 5   2.   **Akane A, Matsubara H, Nakamura S, Takahashi S, Kimura K.** 1994. Identification of  
6       the heme compound copurified with deoxyribonucleic acid (DNA) from bloodstains, a  
7       major inhibitor of polymerase chain reaction (PCR) amplification. *J Forensic Sci.* **39**:  
8       362-372.
- 9   3.   **Chen J-H, Lu F, Lim CS, Kim J-Y, Ahn H-J, Suh I-B, Takeo S, Tsuboi T,**  
10       **Sattabongkot J, Han E-T.** 2010. Detection of *Plasmodium vivax* infection in the Republic  
11       of Korea by loop-mediated isothermal amplification (LAMP). *Acta Trop.* **113**: 61-65.
- 12   4.   **Cox-Singh J, Singh B.** 2008. Knowlesi malaria: newly emergent and of public health  
13       importance. *Trend Parasitol.* **24**: 406-410.
- 14   5.   **Cox-Singh J, Davis TME, Kim-Sung L, Shamsul SSG, Matusop A, Ratnam S,**  
15       **Rahman HA, Conway DJ, Singh B.** 2008. *Plasmodium knowlesi* malaria in humans is  
16       widely distributed and potentially life-threatening. *Clin Infect Dis.* **46**: 165-171.

- 1    6.   **Grab DJ, Lonsdale-Eccles J, Inoue N.** 2005. LAMP for tedpoles. *Nature Methods* **2**:635
- 2            **Han ET, Watanabe R, Sattabongkot J, Khuntirat B, Sirichaisinthop J, Iriko H, Jin L,**
- 3            **Takeo S, Tsuboi T.** 2007. Detection of four *Plasmodium* species by genus- and species-
- 4            specific loop-mediated isothermal amplification for clinical diagnosis. *J Clin Microbiol.*
- 5            **45**:2521-2528.
- 6    7.   **Han ET, Watanabe R, Sattabongkot J, Khuntirat B, Sirichaisinthop J, Iriko H, Jin L,**
- 7            **Takeo S, Tsuboi T.** 2007. Detection of four *Plasmodium* species by genus- and species-
- 8            specific loop-mediated isothermal amplification for clinical diagnosis. *J Clin Microbiol.*
- 9            **45**:2521-2528.
- 10   8.   **Karanis P, Ongerth J.** 2009. LAMP- a powerful and flexible tool for monitoring
- 11            microbial pathogens. *Trend Parasitol.* **25**:498-499.
- 12   9.   **Kawai S, Hirai M, Haruki K, Tanabe K, Chigusa Y.** 2009. Cross-reactivity in rapid
- 13            diagnostic tests between human malaria and zoonotic simian malaria parasite *Plasmodium*
- 14            *knowlesi* infections. *Parasitol Int.* **58**; 300-302.
- 15   10. **Mori Y, Nagamine K, Tomita N, Notomi T.** 2001. Detection of loop-mediated isothermal
- 16            amplification reaction by turbidity derived from magnesium pyrophosphate formation.
- 17            *Biochem Biophys Res Commun.* **289**:150-154.

- 1     11. **Mori Y, Notomi T.** 2009. Loop-mediated isothermal amplification (LAMP): a rapid  
2       accurate, and cost-effective diagnostic method for infectious diseases. *J Infect Chemother.*  
3       **15:** 62-69.
- 4     12. **Nakagami K, Hase T, Notomi T.** 2002. Accelerated reaction by loop-mediated isothermal  
5       amplification using loop primers. *Mol Cell Probes.* **16:** 223-229.
- 6     13. **Notomi T, Okayama H, Masubuchi H, Yonekawa T, Watanabe K, Amino N, Hase T.**  
7       2000. Loop-mediated isothermal amplification of DNA. *Nucleic Acids Res.* **28:** e63
- 8     14. **Poon LLM, Wong BWY, Ma EHT, Chan KH, Chow LMC, Abeyewickreme W,**  
9       **Tangpukdee N, Yuen KY, Guan Y, Looareesuwan S, Peiris M.** 2006. Sensitive and  
10      inexpensive molecular test for falciparum malaria: detecting *Plasmodium*  
11      *falciparum* DNA directly from heat-treated blood by loop-mediated isothermal  
12      amplification. *Clin Chem* **52:** 303-306.
- 13    15. **Singh B, Bobogare A, Cox-Singh J, Snounou G, Adullah HS, Rahman HA.** 1999. A  
14      genus- and species-specific nested polymerase chain reaction malaria detection assay for  
15      epidemiologic studies. *Am J Trop Med Hyg.* **60:** 687-692.
- 16    16. **Singh B, Kim-Sung L, Matusop A, Radhakrishnan A, Shamsul SSG, Cox-Singh J,**  
17      **5Thomas A, Conway DJ.** 2004. A large focus of naturally acquired *Plasmodium knowlesi*  
18      infections in human being. *Lancet.* **363:**1017-1024.

- 1 17. **White NJ.** 2008. *Plasmodium knowlesi*: The fifth human malaria parasite. Clin Infect Dis.
- 2 46: 172-173.

## 1 **Figure legends**

2 Figure 1: Locations and sequences of LAMP targets and priming sites for *P. knowlesi*  $\beta$ -tubulin  
 3 gene. (A) Locations of priming sites of *Pk*LAMP primer set in the reference sequence  
 4 (GenBank accession number: AY639984) are indicated by arrows. (B) Primer sets used for  
 5 amplification of *P. knowlesi*  $\beta$ -tubulin gene in LAMP.

6  
 7 Figure 2: Specificity of *Pk*LAMP for *P. knowlesi*. Panel A, Agarose gel electrophoresis of  
 8 LAMP products from genomic DNA of 13 *Plasmodium* spp. and ethidium bromide staining.  
 9 Panel B, visual detection of LAMP products under UV light using Loopamp fluorescent  
 10 detection reagent. Lane M indicates 200-bp ladder size markers; lane 1, *P. knowlesi*; lane 2, *P.*  
 11 *falciparum*; lane 3, *P. malariae*; lane 4, *P. vivax*; lane 5, *P. ovale*; lane 6, *P. coatneyi*; lane 7, *P.*  
 12 *cynomolgi*; lane 8, *P. inui*; lane 9, *P. simiovale*; lane 10, *P. fieldi*; lane 11, *P. fragile*; lane 12, *P.*  
 13 *gonderi*; lane 13, *P. hylobati*.

14  
 15 Figure 3: Comparison of sensitivities among 3 methods of detection of *Pk*LAMP and  
 16 conventional single-PCR for the detection of *P. knowlesi*  $\beta$ -tubulin gene. Template DNA was  
 17 prepared on serial dilutions of plasmid DNA ( $10^8$  to 1 copies per reaction) containing a  
 18  $\beta$ -tubulin gene for each assay. Panel A, Real-time LAMP assay was monitored by real-time

1 measurement of turbidity. Panel B, Agarose gel electrophoresis of LAMP products. Panel C,  
2 visual detection of LAMP products under UV light using Loopamp fluorescent detection  
3 reagent. Panel D, Agarose gel electrophoresis of single-PCR products using F3 and B3 primers.  
4 Lane M indicates 200-bp ladder size markers (Panel A) and 100-bp ladder size markers (Panel  
5 B); lane 1 to 9,  $10^8$  to 1 copies of plasmid; lane 10, distilled water (Panel B-D).  
6  
7 Table 1: Comparison *Pk*LAMP and nPCR and microscopic examination for detection of *P.*  
8 *knowlesi* using two infected monkeys  
9


**Table 1**

day after infection	Monkey J58 ( <i>P.k.</i> , H strain)					Monkey J64 ( <i>P.k.</i> , Hackeri strain)				
	Parasitemia (%)	<i>Pk</i> LAMP (DNA ex.)	<i>Pk</i> LAMP (w-blood)	nPCR (DNA ex.)	nPCR (w-blood)	Parasitemia (%)	<i>Pk</i> LAMP (DNA ex.)	<i>Pk</i> LAMP (w-blood)	nPCR (DNA ex.)	nPCR (w-blood)
day 0	-	-	-	-	-	-	-	-	-	-
day 1	<0.01	+	+	+	-	-	+	+	-	-
day 2	0.2	+	+	+	-	-	+	+	-	-
day 3	10.8 (autopsy)	+	+	+	+	-	+	+	+	-
day 4						-	+	+	+	-
day 5						-	+	+	+	-
day 6						0.01	+	+	+	-
day 7						0.1	+	+	+	-
day 8						2.0	+	+	+	-
day 9						58.0 (autopsy)	+	+	+	+

DNA ex.: DNA extract      w-blood: whole blood



Figure 1

**B**

Primer	Type	Sequence *
Knowlesi-F3	F3	5' –GGGTGTACACGAACAAGTGT–3'
Knowlesi-B3	B3c	5' –TATCGAGCGTGCGAAGGA–3'
Knowlesi-FIP	F1c-F2	5' –GGTACTCCACCTGTGTATACGCACAGAA <b>TT</b> CGCCTGCATGTATGTATACCAG–3'
Knowlesi-BIP	B1-B2c	5' –GTGAGTGCAACAATGTGAGGAGGT <b>GAATT</b> CCGGAATAAGCATAGGTGGAACC–3'
Knowlesi-FLP	Loop F	5' –GCGCAGACTAGTGCAAGCA–3'
Knowlesi-BLP	Loop B	5' –ACTTCAGTATCCCCCATGTGA–3'

\*Underlining indicates a restriction enzyme site of *Eco R I*.

Figure 2




Figure 3

